




# Disposable Hot Biopsy Forceps



## Reliability at its best

The Disposable Hot Biopsy Forceps made by Micro-Tech allow a reliable and easy procedure to take significant biopsies by high frequency currents. The Hot Disposable Biopsy forceps meet the appropriate IEC Electrical Safety Standards.

- High-class rivet-hinge design
- Excellent price-performance-ratio
- Ergonomic handle
- Good cutting feature and perfect handling



## Coding Rule of IOC

### Note

Note: IOC stands for International Ordering Code.

**HB F XX - X X X XX XXX**  
**(1) (2) (3) (4) (5) (6)**

### Explanation:

- (1) Part **HB F XX**/(1) stands for Hot Biopsy Forceps  
HBF05: oval jaws Hot Biopsy Forceps, with coating
- (2) Part **X**/(2) stands for the size of jaws opening  
1.8 Hot Biopsy Forceps- 1: 4.5mm  
2.3 Hot Biopsy Forceps- 0: 5.5mm; 1: 6.7mm
- (3) Part **X**/(3) stands for the structure type  
0: twin wire 1: four-bar linkage
- (4) Part **X**/(4) stands for the spike  
0: without spike
- (5) Part **XX**/(5) stands for the diameter of jaws  
18: 1.8mm 23: 2.3mm
- (6) Part **XXX**/(6) stands for the working length  
120: 1200mm 160: 1600mm 180: 1800mm 230: 2300mm

### For example:

HBF05-11023230 refers to oval jaws Hot Biopsy Forceps with coating, the size of jaws opening is 6.7mm, 4-bar linkage, without spike, the diameter of jaws is 2.3mm, and the working length is 2300mm.

## Specification

IOC*	EOC*	Dia. (mm)	Working Length (mm)	Spike	Coating	Jaws	Jaws Opening (mm)
HBF05-11023160	-----*	2.3	1600	No	Yes	oval	6.7
HBF05-01023160	-----	2.3	1600	No	Yes	oval	5.5
HBF05-11023180	-----	2.3	1800	No	Yes	oval	6.7
HBF05-01023180	-----	2.3	1800	No	Yes	oval	5.5
HBF05-11023230	-----	2.3	2300	No	Yes	oval	6.7
HBF05-01023230	-----	2.3	2300	No	Yes	oval	5.5
HBF05-11018120	-----	1.8	1200	No	Yes	oval	4.5
HBF05-11018160	-----	1.8	1600	No	Yes	oval	4.5
HBF05-11018180	-----	1.8	1800	No	Yes	oval	4.5

.....

### Note

IOC\*: Stand for International Ordering Code.

EOC\*: Stand for European Ordering Code.

-----\*: Stand for EOC is the same as IOC.

