



Rectal Stent



The perfect solution for the rectum

The Rectal Stent is preloaded on OTW (over-the-wire) delivery system. The specific end design has the perfect fit for the anatomic and peristaltic conditions.

- Made of single Nitinol wire with the "thermal memory effect" allowing the stent to expand to its preset shape after released
- Mesh structure and atraumatic ends to reduce irritations to the patient
- Strong radial force and anatomy-fit design to avoid migration
- Extraction threads for stent repositioning and/or removal
- Radiopaque marks for easy orientation and exact implantation under X-ray
- Custom-made products available



Specification

IOC*	EOC*	Body Dia. (mm)	Flange Dia. (mm)	Length (mm)	Covering (mm)
Stents without covering					
ST44-331.25.090	ST04-111-25.090	25	30	090	Without covering
ST44-331.30.100	ST04-111-30.100	30	36	100	Without covering
ST44-331.23.120	ST04-111-30.120	30	36	120	Without covering
Stents with partial covering					
ST44-332.25.090	ST04-112-25.090	25	30	090	70
ST44-332.30.100	ST04-112-30.100	30	36	100	80
ST44-332.30.120	ST04-112-30.120	30	36	120	100
Stents with end-to-end covering					
ST44-334.25.090	ST04-113-25.090	25	30	090	090
ST44-334.30.100	ST04-113-30.100	30	36	100	100
ST44-334.30.120	ST04-113-30.120	30	36	120	120

.....

For more specifications, please refer to CE range list of intestinal stent. (See Page 55).


Note

IOC*: Stand for International Ordering Code.


EOC*: Stand for European Ordering Code.

Delivery System


Ref	Dia. (mm)	Length (mm)	Guide Wire
OTW (Over the wire) System	8, 8.5	700	0.035"




Olive tip



Front handle



Loop assembly



Balloon handle